



Production line



Sterilization chamber



Clean Room

# FMI Medical



## Utility face mask

Produced, Sterilized in our clean room and packaged in the USA

With the help of FMI engineers, designers and buyers, FMI has converted part of our facility headquartered in Solon, Ohio to produce Utility Flat and KN95 face masks to help keep those on the front lines, their families and the people they service safe.

### Product facts:

- Single use face masks
- 3 layers combination of nonwoven and melt blown/spray fabric
- Poly coated nose wire
- Elastic ear Loops
- Made in Cleveland, Ohio

### Material facts:

- Outside layer: 25gsm filament nonwoven fabric
- Middle layer: 25gsm melt spray fabric
- Inner layer: 25gsm hygroscopic nonwoven fabric
- Nose band: Poly coated metal material
- Ear loop: Spandex


\*The specifications and descriptions above are based on the data provided by our raw material suppliers.

### Intended use:

A physical barrier for general public daily use.

### Warning:

Materials in direct contact with the skin may cause allergic reactions in some sensitive individuals.



**Manufacturers:** FMI Medical Systems, Inc.  
**Business address:** 29001 Solon Rd., Unit A, Solon, OH 44139, USA  
**Telephone number:** 877-757-1195  
**Email:** FMIMED.Marketing@fmimedical.com  
**Website:** www.fmimedical.com

# FMI Medical

## Utility face mask

### Caution:

Store at room temperature in cool and dry place. Do not use if the packing is damaged. Manufacture date/lot and expiration date, see packing box or label inside.

### FDA COVID-19 EUA Guidance:

- This face mask is NOT a surgical mask. It is NOT to provide liquid barrier protection.
- This face mask is NOT for use in a clinical setting where the infection risk level through inhalation exposure is high.
- This mask is NOT for antimicrobial or antiviral protection or related uses for infection prevention or reduction or related uses.
- This face mask is NOT a respiratory protective device.
- This face mask is NOT for high risk aerosol-generating procedures.
- Shall not be used in the presence of a high intensity heat source or flammable gas.

### Certifications:

FMI Medical is working on FDA EUA and long term 510k.

### Product options:

Part Number	Description	Color	Layers
40-0001	Utility face mask, 10 pcs package	Blue	3
40-0005	Utility face mask, 10 pcs package	Pink	3
40-0006	Utility face mask, 50 pcs package	Blue	3
40-0007	Utility face mask, 50 pcs package	Pink	3

### Products photos:

- 10 pcs package is plastic bag
- 50 pcs package is paper carton
- Offering shipping box for 1200 pcs order



**Manufacturers:** FMI Medical Systems, Inc.  
**Business address:** 29001 Solon Rd., Unit A, Solon, OH 44139, USA  
**Telephone number:** 877-757-1195  
**Email:** FMIMED.Marketing@fmimedical.com  
**Website:** www.fmimedical.com